

LIQUID BASED MICROBIOLOGY

FecalSwab™

COLLECTION & PRESERVATION
SYSTEM FOR ENTERIC DIAGNOSTICS



Innovating Together™

copangroup.com




FecalSwab™ is a **Collection, Preservation and Transport** system for Feces and Rectal samples.

It is specifically designed for **enteric pathogens**.



KEY FEATURES

- › Ready-to-use
- › Sample **preservation**: 48h at RT and 72h at 4°C**
- › **Compact container, neat, easy** to transport, use and store
- › **Versatile**: sample collected from both defecated stool and rectal swab
- › **Validated** for: Culture
- › **Supported** analyses: Microscopy *, rapid antigen*, PCR*
- › Tested with **molecular platforms***** : Biofire Film Array, BD Max, R Bio-Pharm, GenomEra
- › Tube with **maximum sample quantity** indication
- › FLOQSwab™ with molded **breaking point**
- › Tube with **capture cap**
- › **New package**: Vi-pack



* Analyses applicable after user 's validation


** For *C. difficile* 24h at RT and 48h at 4°C

*** Publications available upon request.


SAMPLE WORKFLOW

Microscopy°

Bacteria



Parasites



Culture

E.coli
Escherichia coli O157:H
Shigella flexneri
S. typhimurium
E.faecalis (VRE)

*V.para*haemoliticus
Y. enterocolitica
C. jejuni
C. difficile



Tested at RT and 4°C

Manual & Automation


Antigen detection°

Adenovirus - *Rotavirus* - *C. difficile*



Molecular assay°

Nucleic acid extraction° & PCR°



°Analyses applicable after user's validation. Publications available upon request.



FecalSwab™

ORDERING INFORMATION

Code	Product Description	Pack Size
470CE	FecalSwab™ for manual use, 12x80mm tube filled with 2 ml Modified Cary Blair medium + 1 regular FLOQSwab™	300 pieces per carton (6 boxes of 50 pieces)
470CE.A	FecalSwab™ for automation, 12x80mm tube filled with 2 ml Modified Cary Blair medium + 1 regular FLOQSwab™	300 pieces per carton (6 boxes of 50 pieces)
4E020S.A	FecalSwab™ for automation, 12x80mm rounded bottom tube, filled with 2 ml Modified Cary Blair medium + 1 regular FLOQSwab™	300 pieces per carton (6 boxes of 50 pieces)
4C028S*	FecalSwab™ for manual use, 12x80mm tube filled with 2 ml Modified Cary Blair medium + 1 rectal FLOQSwab™	300 pieces per carton (6 boxes of 50 pieces)
4U031S	FecalSwab™ in bulk for stool specimen: FLOQSwab™ stool transfer device in resealable pack + 12x80mm tube filled with 2 ml Modified Cary Blair medium in bulk	300 pieces per carton (6 boxes with 50 pieces of tubes + 1 resealable pack of 50 stool transfer device)

*under development



- › European Patent #1608268
- › US Patent #8,114,027
- › US Patent #8,317,728
- › Canadian Patent #2,515,205
- › Japanese Patent #4579902
- › Chinese Patent #101103931
- › Australian Patent #2004226798
- › New Zeland Patent #541560

This document contains information on products which is targeted to a wide range of audiences and may contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin.



For general distribution enquiries:

COPAN ITALIA SPA
 via F. Perotti 10
 25125 - Brescia, Italy
 t. +39 030 2687211
 info@copangroup.com
 www.copangroup.com

For diagnostic companies enquiries:

COPAN FLOCK TECHNOLOGIES SRL
 c/o Futura Science Park
 via F. Perotti 18
 25125 - Brescia, Italy
 t. +39 030 3666100
 info@copangroup.com
 www.copangroup.com

Brochure hardcopies may not include the latest updates and changes. Please refer to Copan website to view and/or download the most recent controlled version of the brochure. Brochures are mainly intended for marketing purposes. Always consult product inserts, instructions for use and manuals for the appropriate use of the products. ©2016 Copan Italia. All rights reserved. The trademarks mentioned herein are property of Copan.